

Fracture Implants



MECHANICALLY STABLE

STAINLESS STEEL

POLYAXIAL HOLE



Liberty Lock Trauma™ plates can be pre-contoured in all planes yet still have rigidity close to limited contact plates. Between each screw location are bending zones to control where the bend occurs so not to damage the locking holes to ensure proper fixation.

▶ Surgical Options

Our plates are available in 2.4, 2.7, and 3.5mm sizes.

▶ Flexibility

This innovative hole design allows for polyaxial fixation. Place and lock screws where desired.

LIBERTY LOCK TRAUMA TPLO PRODUCT BENEFITS:

- ▶ Compatible with industry standard instruments and implants
- ▶ Multi-plane contouring
- ▶ 6 degrees of freedom
- ▶ Manufactured at ISO certified facilities in the UK and Germany

Ask about our complete line of instrumentation and sterile procedural packs.